

SIGMA F.O. LED OTOSCOPE 2.5V with rechargeable handle and battery - pouch - black



Product Code	31580
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2831178
CND	Z12149006
EAN/UPC	8023279315806
GTIN	08023279315806
GMDN	12849

Description

- SIGMA F.O. LED OTOSCOPE 2.5V - rechargeable handle with battery
- Black

Rechargeable handle can be used with "AA" standard batteries or rechargeable battery code 31587 (needs charging station code 31586)

Full battery charge in 2.5 hours. Charging station can be connected to PC or laptop via standard Type C USB cable (not included)

High quality pocket diagnostic instruments designed for professional use in ENT and ophthalmoscopy. Light weight and manufactured from impact resistant re-enforced plastic.

Manufactured with high quality materials and top-notch workmanship ensures excellent reliability in every examination.

Battery life with Led bulb 6 time longer than with traditional bulb

- distal fibre optic illumination with super bright LED light
- swivel type window with 3X magnification lens for sharp images and minimal reflection
- insufflation port for pneumatic testing of tympanic mobility

Connector port (code 31532) and insufflation bulb (code 31500) should be ordered separately.