

SIGMA F.O. LED OTO-OPHTHALMOSCOPE SET with 2 rechargeable handles - case



Product Code	31584
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2734775
CND	L1699
EAN/UPC	8023279315844
GTIN	08023279315844
GMDN	45084

Description

Oto-ophthalmoscope set with 2 rechargeable handles and batteries

Type: LED
Packing: Case

This set is supplied with oto and ophthalmo head, 2 handles, 14 disposable specula (7 with Ø 2.5 mm, 7 with Ø 4 mm), 2 rechargeable batteries.