

SIGMA-C LED OTOSCOPE - black - pouch



Product Code	31591
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2831178
CND	Z12149006
EAN/UPC	8023279315912
GTIN	08023279315912
GMDN	12849

Description

SIGMA-C LED OTOSCOPE - black - pouch

- made of light weight impact resistant re-enforced plastic and highly corrosion resistant stainless steel material for exterior finish
- swivel type window with 3X magnification lens for sharp image and minimal reflection
- optional insufflation port for pneumatic testing of tympanic mobility (needs codes 31532 and 31483)
- compatible with Gima ear specula 31485-8, 31491-2, 31495-9
- 50.000 hours lifetime
- super bright 25,000 LUX
- extended battery life of 50 hours
- colour temperature 5,800 K, CRI>78 ensures real colour.

Supplied in a pouch with 14 specula (7x2.5 mm, 7x4.0 mm) in transparent tube

Multilingual manual: GB, FR, IT, ES, PT, DE, GR, RO, PL, Arabic.