

VISIO 5GE HQ BASIC OTO-OPHTHALMO WALL DIAGNOSTIC STATION - LED



Product Code	31621
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2678491
CND	V9099
EAN/UPC	8961102783241

Description

Visio 5GE HQ Basic oto-ophthalmo diagnostic station - LED

Include power transformer, 2 handles, oto-ophthalmo heads, dispenser with 100 disposable grey speculum (50 Ø 2.5 mm, 50 Ø 4 mm), wall board for easy mounting or removal of diagnostic station.

Modern design diagnostic stations to satisfy all practitioners and clinics needs.

Visio 5GE HQ Basic oto-ophthalmo diagnostic station - LED
31622 Visio 5XE HQ oto-ophthalmo diagnostic station - Xenon

Include power transformer, 2 handles, oto-ophthalmo heads, dispenser with 100 disposable grey speculum (50 Ø 2.5 mm, 50 Ø 4 mm), wall board for easy mounting or removal of diagnostic station.

Modern design diagnostic stations to satisfy all practitioners and clinics needs.

Regulated output prolongs the bulb life and guarantees optimal light intensity

- ergonomic designed non-slip handles
- constant output transformer with integrated surge protection
- soft start technology to introduce voltage gently, eliminates cold-shock when turning ON, increasing bulb life
- adjustable light intensity, low power consumption
- automatic ON-OFF switch once the handle is removed or replaced in the cradle
- precise brightness control, between 3 and 100%
- spiral cord of 3 meters

Supplied with multilingual manual: GB, FR, IT, ES.

Technical Specifications

Apertures: 7
Focusing lenses: 28
Diopters range: -25 to +40