

# HEINE MINI 3000 LED DIRECT OTOSCOPE with reusable tips - blue - D-008.70.220



Product Code	31663
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

EAN/UPC	4053755202607
---------	---------------

## Description

### HEINE MINI 3000 LED 2.5V OTOSCOPE

Mini 3000 LEDHQ direct illumination otoscope

- Blue
- Code HEINE: D-008.70.220
- Reusable Tips 4 pcs x 4 mm
- AllSpec disposable Tips 10 pcs x 4 mm
- Battery handle (Batteries included)

Compact otoscope with direct LED illumination. Natural colour rendering, no need to ever change the LED. Exclusive design.

#### Detailed view

- it eliminates reflections due to the head's black inner coating
- high-resolution images
- colour temperature 4,000°K
- 3x magnification
- particularly bright due to LEDHQ

#### Intuitive operability and practical attachment clip

- it switches off automatically as soon as it is clipped to the pocket
- easy instrumentation. Viewing window fixed, can be swivelled to the side
- operating time up to 10 hours