

HEINE MINI 3000 F.O. LED OTOSCOPE with reusable tips - black - D-008.70.120



Product Code	31674
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

EAN/UPC 4053755158676

Description

HEINE MINI 3000 LED FIBER OPTICS (F.O.) OTOSCOPES

Mini 3000 LEDHQ F.O. otoscope

- Black
- Code HEINE: D-008.70.120
- Reusable Tips 4 pcs x 4 mm
- AllSpec disposable Tips 10 pcs x 4 mm
- Battery handle (Batteries included)

Compact otoscope with especially bright illumination and authentic colour rendering due to LEDHQ (no need to ever change the LED). Exclusive design and sturdy housing.

Detailed view

- high-resolution images, minimal reflection
- 3x magnification
- reflex-free illumination of tympanum and ear canal due to fiber optic technology
- particularly bright due to LEDHQ
- operating time up to 10 hours

Intuitive operability and practical attachment clip

- it switches off automatically as soon as it is clipped to the pocket
- easy instrumentation. Viewing window fixed, can be swivelled to the side
- pneumatic eardrum testing by means of an integrated insufflation port