

HEINE MINI 3000 F.O. LED OTOSCOPE with reusable tips - black - D-008.70.120



Product Code	31674
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

EAN/UPC	4053755158676
---------	---------------

Description

HEINE MINI 3000 LED FIBER OPTICS (F.O.) OTOSCOPES

- Mini 3000 LEDHQ F.O. otoscope
- Black
 - Code HEINE: D-008.70.120
 - Reusable Tips 4 pcs x 4 mm
 - AllSpec disposable Tips 10 pcs x 4 mm
 - Battery handle (Batteries included)

Compact otoscope with especially bright illumination and authentic colour rendering due to LEDHQ (no need to ever change the LED). Exclusive design and sturdy housing.

- Detailed view
- high-resolution images, minimal reflection
 - 3x magnification
 - reflex-free illumination of tympanum and ear canal due to fiber optic technology
 - particularly bright due to LEDHQ
 - operating time up to 10 hours

- Intuitive operability and practical attachment clip
- it switches off automatically as soon as it is clipped to the pocket
 - easy instrumentation. Viewing window fixed, can be swivelled to the side
 - pneumatic eardrum testing by means of an integrated insufflation port