



|               |                        |
|---------------|------------------------|
| Product Code  | 31677                  |
| Unit of sale  | 1 pc                   |
| Minimum order | 1                      |
| Type          | Kit of medical devices |
| Class         |                        |

|         |               |
|---------|---------------|
| UK-REP  | Yes           |
| EAN/UPC | 4053755158645 |

## Description

### HEINE MINI 3000 LED FIBER OPTICS (F.O.) OTOSCOPES

Mini 3000 LEDHQ F.O. otoscope

- Blue
- Code HEINE: D-008.70.110
- AllSpec disposable Tips 10 pcs x 4 mm
- Battery handle (Batteries included)

Compact otoscope with especially bright illumination and authentic colour rendering due to LEDHQ (no need to ever change the LED). Exclusive design and sturdy housing.

#### Detailed view

- high-resolution images, minimal reflection
- 3x magnification
- reflex-free illumination of tympanum and ear canal due to fiber optic technology
- particularly bright due to LEDHQ
- operating time up to 10 hours

#### Intuitive operability and practical attachment clip

- it switches off automatically as soon as it is clipped to the pocket
- easy instrumentation. Viewing window fixed, can be swivelled to the side
- pneumatic eardrum testing by means of an integrated insufflation port