



Product Code	31688
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	Yes
EAN/UPC	4053755158737

## Description

### HEINE MINI 3000 LED 2.5V OPHTHALMOSCOPES

#### Mini 3000 LEDHQ ophthalmoscope

- Black
- Code HEINE: D-008.71.120
- Battery handle (Batteries included)

Modern pocket ophthalmoscope with LEDHQ. Concentrated, bright light for optimum illumination, a brilliant image and authentic colour rendering. Fixation star aperture with grey filter to reduce light intensity. Excellent optics and outstanding colour rendering with the smallest dimensions.

#### Reliable diagnosis

- simple examination process possible without pupil dilation
- clear, high-resolution image with minimal reflections
- 18 lenses ranging from - 20 to + 20 D prevent higher refractive errors
- fixation star including grey filter for patients who are sensitive to light
- 5 different apertures

#### Intuitive operability

- it switches off automatically as soon as it is clipped to the pocket
- practical one-finger operation
- handle and head are in 2 parts for easy servicing and expandability