

# HEINE MINI 3000 F.O. OPHTHALMOSCOPE - rechargeable handle - black



Product Code	31694
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	2779683
CND	Z12120114
EAN/UPC	8023279316940

## Description

### MINI 3000 OPHTHALMOSCOPES 2.5 V - RECHARGEABLE HANDLE

- F.O. Xenon
- RECHARGEABLE HANDLE (Need Heine Mini 3000 NT charger code 31692)
- Colour: Black
- HEINE code: D-001.71.200

Modern pocket ophthalmoscope with fixation star aperture. Excellent optics and best colour rendering in minimum size.

- XHL Xenon Halogen Technology. Concentrated bright light for perfect illumination, a brilliant image and accurate colour rendering
- the HEINE optics produce crisp, clear images of the whole area under examination and minimize the reflections
- 5 different apertures. A full-function instrument incl. fixation star
- grey filter for light sensitive patients
- range of lenses: 18 lenses from -20 D to +20 D
- easy workflow, because ophthalmoscopy is possible without dilation of the pupils
- attachment clip with integrated on/off switch
- replaceable batteries. Size AA or optional rechargeable battery (31693) with mini NT table charger (31692)