

# HEINE MINI 3000 LED F.O. OTO + OPHTHALMO SET - black - D-886.11.021



Product Code	31696
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	Yes
EAN/UPC	4053755161256

## Description

### HEINE MINI 3000 LED DIAGNOSTIC SET - OPHTHALMOSCOPIES AND OTOSCOPIES

Mini 3000 LEDHQ F.O. oto-ophthalmoscope

- Black
- Code HEINE: D-886.11.021
- Reusable Tips 4 pcs x 4 mm
- AllSpec disposable Tips 10 pcs x 4 mm
- Battery handle (Batteries included)
- Hard case

Modern pocket-size ophthalmoscope in combination with a compact and very bright otoscope. Handy and with LEDHQ illumination. Concentrated, bright light for optimum illumination, a brilliant image and authentic colour rendering.