

HEINE MINI 3000 LED DIRECT OTO + OPHTHALMO SET - black -D-872.11.021



Product Code	31698
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

EAN/UPC	4053755202560
---------	---------------

Description

HEINE MINI 3000 LED DIAGNOSTIC SET - OPHTHALMOSCOPES AND OTOSCOPES

Mini 3000 LEDHQ direct illum. oto-ophthalmoscope

- Black
- Code HEINE: D-872.11.021
- Reusable Tips 4 pcs x 4 mm
- AllSpec disposable Tips 10 pcs x 4 mm
- Battery handle (Batteries included)
- Hard case

Modern pocket-size ophthalmoscope in combination with a compact and very bright otoscope. Handy and with LEDHQ illumination. Concentrated, bright light for optimum illumination, a brilliant image and authentic colour rendering.