

# HEINE BETA 200S LED OPHTHALMOSCOPE - 3.5V C-261.28.388



Product Code	31737
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	1725637
CND	Z12019080
EAN/UPC	4053755184330
UMDNS	12-817

## Description

HEINE BETA 200S LED OPHTHALMO - 3.5 V with USB rechargeable handle - C-261.28.388

Superior aspherical optics and 74 single-diopter steps

- optical system: corneal and iris reflexes are eliminated. A full view is guaranteed even with small pupils.
- XHL (Xenon Halogen Technology) or LED: bright, white light.
- recessed, multi-coated viewing window: avoids stray light.
- optics mounted on precision metal chassis.
- ergonomic shape: instrument fits the orbita comfortably in any position.
- 7 apertures with separate red-free filter:
- lens range: perfect focus even with high refractive errors.
- exact colour rendering

Made in Germany.