

HEINE BETA 200 F.O. LED OPHTHALMOSCOPE HEAD -
2.5 V - C-001.30.100



Product Code	31743
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

CND	Z12120114
EAN/UPC	4053755119981

Description

HEINE code: C-001.30.100
OPHTHALMOSCOPES HEADS: BETA 200 LED ophthalmoscope head - 2.5 V
Compatibility: Battery handles