

HEINE BETA 200 F.O. LED OPHTHALMOSCOPE HEAD - 3.5 V - C-008.30.100



Product Code	31744
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

CND	Z12120114
EAN/UPC	4053755182534

Description

HEINE code: C-008.30.100
OPHTHALMOSCOPES HEADS: BETA 200 LED ophthalmoscope head - 3.5 V
Compatibility: Rechargeable handles