



Product Code	31745
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

RDM (NSIS)	1717919
CND	Z12120118
EAN/UPC	4053755121533
UMDNS	13-372
PARAF	965528425

Description

HEINE BETA 200 RETINOSCOPE - slit - 2.5 V - C-034.10.118

With ParaStop® for finding the neutralization point (US pat. 5,859,687): precise, easy selection of a parallel beam

- dustproof: maintenance free
- single control for vergence and rotation: comfortable operation
- Integrated polarisation filter: eliminates stray light and internal reflexes for a brighter pupilling reflex
- streak or spot retinoscope: simply by changing the bulb
- XHL (xenon halogen technology): bright, white light. Bright fundus reflex, easy recognition of neutralization point
- ergonomic shape: protects the examiner's orbit from stray light

Supplied in hard case with one spare bulb.