

# HEINE BETA 200 LED RETINOSCOPE 2.5V - slit - C-034.10.118



Product Code	31746
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC	4053755121533
PARAF	965576162

## Description

HEINE BETA 200 LED RETINOSCOPE - slit - 2.5 V - C-034.10.118

Hand-held retinoscope for objective refraction with totally even LEDHQ illumination for a precise streak. It enables quick and easy recognition of the fundus reflex.

Reliable diagnosis: precise adjustment of the parallel beam path due to ParaStop. The LEDHQ provides a particularly bright, 1.1 mm wide streak and a clear fundus reflex, with patented 3-100% dimming and a polarisation filter to eliminate glare and reflections.

Enhanced durability: maintenance-free, dustproof design with robust metal controls. Includes a virtually unlimited LED lifespan.

Ergonomic operation: intuitive one-finger control for both vergence and line rotation.

Supplied with battery handle, 1 spare bulb and hard case.

## Technical Specifications

Colour temp.: 3000 K

Length of the streak: 35 mm at 500 mm distance

Width of the streak: 1.1 mm at 500 mm distance

Working distance: 500 mm

Size - weight: 129x37x43 mm - 90 g