

# HEINE MPACK UNPLUGGED RECHARGEABLE BATTERY for Omega 500 and ML4 Headlight - X-007.99.665



Product Code	31752
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

EAN/UPC	4053755144723
---------	---------------

## Description

RECHARGEABLE BATTERIES for omega 500 and ml4 headlight

- HEINE code: X-007.99.665
- mPack unplugged rechargeable battery