

HEINE 049 BULB 3.5V for Alpha+, Beta100, K100, K100 Vet otoscopes



Product Code	31787
Unit of sale	1 pc
Minimum order	1
Type	Spare part
Class	

UK-REP	Yes
--------	-----

EAN/UPC	4053755143849
---------	---------------

Description

HEINE bulb Xenon 3.5 V

Heine ref.: 049

Compatible with following HEINE products: F.O. Alpha+, Beta 100, K 100, Beta 100 vet otoscopes; Alpha+, Beta® tongue holders