

RI-FORMER XENON DIAGNOSTIC STATION - 3,5 100-240V - Standard



Product Code	31800
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
EAN/UPC	8023279318005

Description

MODULAR DIAGNOSTIC STATION - XENON - 100-240V, 50/60 Hz
• Ri-Former Diagnostic station 3.5 -standard (include 2 handles, oto and ophthalmo head, ear specula dispenser, clock)

- RI-FORMER MODULAR DIAGNOSTIC STATION
- The modular diagnostic station can be individually combined according to requirements and the room situation, it is variably extendable to the right and the left.
- The basic Ri-former® module is available with one handle and can be supplemented with up to 4 further handle modules.
- 3.5V Xenon or Led illumination creates lighting conditions similar to daylight
 - optimum electronic efficiency using an adjustable dimmer
 - constant output voltage independent from the mains supply prolongs the service life of the lamps and guarantees optimum light intensity
 - automatic switching on/off of handles when removing or replacing them in the handle unit saves energy and increases the service life of the lamps
 - non-slip, ergonomically designed handles with rheostat: infinitely variable regulation of brightness as well as manual switching on/off possible on the handle
 - spiral cord expandable up to 3 m

Made in Germany.