

# RI-FORMER XENON DIAGNOSTIC STATION - 3,5 100-240V - Standard



Product Code	31800
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

EAN/UPC	8023279318005
---------	---------------

## Description

MODULAR DIAGNOSTIC STATION - XENON - 100-240V, 50/60 Hz

• Ri-Former Diagnostic station 3.5 -standard (include 2 handles, oto and ophthamo head, ear specula dispenser, clock)

### RI-FORMER MODULAR DIAGNOSTIC STATION

The modular diagnostic station can be individually combined according to requirements and the room situation, it is variably extendable to the right and the left.

The basic Ri-former® module is available with one handle and can be supplemented with up to 4 further handle modules.

- 3.5V Xenon or Led illumination creates lighting conditions similar to daylight

- optimum electronic efficiency using an adjustable dimmer

- constant output voltage independent from the mains supply prolongs the service life of the lamps and guarantees optimum light intensity

- automatic switching on/off of handles when removing or replacing them in the handle unit saves energy and increases the service life of the lamps

- non-slip, ergonomically designed handles with rheostat: infinitely variable regulation of brightness as well as manual switching on/off possible on the handle

- spiral cord expandable up to 3 m

Made in Germany.