

RIESTER REUSABLE EAR SPECULUM - diam. 4 mm



| | |
|---------------|----------------|
| Product Code | 31819 |
| Unit of sale | box of 10 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|---------|---------------|
| EAN/UPC | 4045396011618 |
|---------|---------------|

Description

AUTOCLAVABLE MINI EAR SPECULUM for all Riester otoscopes

Riester ear speculum Ø 4 mm - black

Ear-speculum are compatible with Ri-scope L1, L2, E-scope, Ri-mini, Pen-scope