

RI-SCOPE L LED OTOSCOPE L2 - 2.5V - AA-handle - 3708



Product Code	31829
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC 4045396140127

Description

RI-SCOPE L LED OTOSCOPE L2 - 2,5 V - AA handle - 3708

Ri-scope®L otoscopes with optical glass and reduced reflection combine excellent optical properties

- optimal beaming and transmission of the light.
- bi-directional swivelling, sealing optical glass 3X
- Suitable for pneumatic test (supplied without ball).

Supplied in deluxe hard case, 4 reusable and 3 disposable ear specula.

Spare ear speculum: original 31817-31819 (autoclavable) or 31827 - 31828 (disposable). Compatible 31485 to 31499 (disposable).

Made in Germany.