


RI-SCOPE L LED OTOSCOPE L2 - 3.5V - C-handle for Ri-Accu L - 3703-550



Product Code	31830
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1298065
CND	Q0399
EAN/UPC	4045396148697
GMDN	36093

Description

RI-SCOPE L LED OTOSCOPE L2 - 3,5 V - C handle - 3703-550

Ri-scope®L otoscopes with optical glass and reduced reflection combine excellent optical properties and unparalleled performance with the reliability and efficiency of LED/lithium-ion technology.

- optimal beaming and transmission of the light.
- bi-directional swivelling, sealing optical glass 3X
- Suitable for pneumatic test (supplied without ball).
- C handle runs Li-ion batteries, not included.

Supplied in deluxe hard case, 4 reusable and 3 disposable ear specula.

Spare ear speculum: original 31817-31819 (autoclavable) or 31827 - 31828 (disposable). Compatible 31485 to 31499 (disposable).

Made in Germany.