



Product Code	31831
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
EAN/UPC	4045396140202

Description

RI-SCOPE L XL OTOSCOPE L2 - 2.5 V - C handle - 3716

- optimal beaming and transmission of the light
- bi-directional swivelling, sealing optical glass 3X
- suitable for pneumatic test (supplied without ball)
- C handle runs Li-ion batteries, not included.

Supplied in deluxe hard case with 10 ear specula.

Include spare bulb for ophthalmoscope only.

Spare ear speculum: original 31817-31819 (autoclavable) or 31827 - 31828 (disposable).

Compatible 31485 to 31499 (disposable)

Made in Germany.