



Product Code	31832
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
EAN/UPC	4045396152991
GMDN	36093

Description

RI-SCOPE® L XL OPHTHALMOSCOPE L2 SET - 3.5 V - C handle - 3725-550

The enhanced basic model with separately engageable redfree, blue and polarization filters for each aperture.

- easy-to-operate aperture hand-wheel with semi-circle, small/medium/large circle, fixation star, and slit
- dioptre wheel with 29 corrective lenses Plus 1-10, 12, 15, 20, 40 Minus 1-10, 15, 20, 25, 30, 35
- optimized high-performance optics with aspherical condenser lens
- parallel beam path

Supplied in deluxe hard case with magnifying glass with 5-fold magnification.
Spare bulb included.

Made in Germany.