

RIESTER BULB 10600 - XL 2.5 V



| | |
|---------------|----------------|
| Product Code | 31854 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

Description

Riester original bulbs.

Riester ref.: 10600

Bulb: XL 2.5V

Compatible with following Riester products: Otoscope ri-mini®, ri-scope® L2/L3 and e-scope®