



Product Code	31863
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC 4045396209855

## Description

RIESTER ELITEVUE LED OTO-ophthalmo 3.5 V - 2210-203

Both with C-handle for rechargeable Li-Ion battery

EliteVue offers twice the field of view and 66% greater magnification (at 15 mm distance) compared to standard otoscopes.

Improved view of auditory canal and tympanic membrane guarantees more efficient diagnosis.

Suitable for throat illumination.

4,000 K LED illumination (CRI: 92).

5X magnification.

Supplied in a case with specula

Made in Germany.