

RIESTER UNI I OTOSCOPE - Vacuum 2.7 V - C-handle - 2010



Product Code	31870
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

EAN/UPC	4045396002739
---------	---------------

Description

UNI® I OTOSCOPE - Vacuum 2.7 V - C-handle - 2010

Durable and user-friendly otoscopes suitable for ENT examination. Equipped with removable magnifying lens with 4-fold magnification that can be swiveled to both sides and removable cover glass for insertion of external instruments.

- sealed system for performing pneumatic otoscopy
- plug-in connector for quick and secure attachment to handle
- ageing-resistant chrome-plated metal head

Supplied with reusable ear specula 2, 3 and 4 mm, battery handle type C in a silver-coloured plastic case with two snap fasteners.