



Product Code	31875
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
EAN/UPC	4045396002760

Description

ECONOM® DIAGNOSTIC SET - Vacuum 2.7 V/XL 2.5 V - C-handle - 2050

Econom® portable diagnostic set for ear, eye, nose and throat treatment (EENT).

It includes Uni® I otoscope head, Uni® II May ophthalmoscope head, nasal speculum expandable, bent arm illuminator, 2 laryngeal mirrors No. 3 and 4, tongue blade holder, replacement lamp otoscope (also for bent arm illuminator).
 Supplied in hard protective deluxe case with sliding safety catches.

Uni® II May-ophthalmoscope with durable ABS plastic head

- plug-in connector for fast and secure attachment to handle
- focusing wheel with corrective lenses from 0 to +20 or 0 to -20 dioptres
- aperture with large circle for examination of the fundus of the eye