



|               |                |
|---------------|----------------|
| Product Code  | 33209          |
| Unit of sale  | box of 500 pcs |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |                |
|------------|----------------|
| UK-REP     | No             |
| CH-REP     | No             |
| RDM (NSIS) | 57506          |
| CND        | C020501        |
| PARAF      | 977813601      |
| GTIN       | 05707480021987 |
| GMDN       | 35035          |

## Description

- AMBU BLUE SENSOR R
- highly conductive wet gel
  - superior adhesion
  - comfortable foam backing
  - large measuring area
  - unique offset connector
  - high quality Ag/AgCl sensor

BlueSensor is the world's most time and cost efficient range of electrodes that delivers unmatched readings, user convenience and patient comfort. BlueSensor ensures that you don't get your workflow disrupted by false alarms or have to waste time resolving problems with electrode adhesion.

Exact monitoring clear picture Blue Sensor electrodes give a fast and reliable trace without the need for skin preparation, maximizing comfort for the patient. We know that fast, exact monitoring that ensures a clear picture of the patients' heartbeat is essential to accurate diagnoses.

## Technical Specifications

- AMBU code: R-00-S/25
- Description: BLUE SENSOR R
- Application: Stress-test
- Gel type: Wet
- Backing material: Foam
- Size(mm): 57 x 48
- Shape: Round with offset
- Shell life\* (month): 24
- Radiolucent: NO
- Connector type: SNAP
- Minimum order: Box of 500 pcs.
- Packaging: 20 pouches of 25 pcs.

\*Shelf life means maximum life from production date. Average Gima shelf life is 50 to 90% of maximum life.