

KLEENSPEC PAEDIATRIC EAR SPECULA 52432-UB - disposable



Product Code	32105
Unit of sale	box of 850 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
UMDNS	13662
GMDN	34897

Description

WELCH-ALLYN KLEENSPEC DISPOSABLE EAR SPECULUM

Colour: black

Ø: 2.75 mm

To be used with: Welch-Allyn, MacroView

Minimum order: box of 850