

# KLEENSPEC ADULT EAR SPECULA 52434-UB - disposable



Product Code	32106
Unit of sale	box of 850 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

EAN/UPC	0732094107654
UMDNS	13662
GMDN	34897

## Description

WELCH-ALLYN KLEENSPEC DISPOSABLE EAR SPECULUM

Colour: black

Ø: 4.25 mm

To be used with: Riester Ri-scope L2, L3

Minimum order: box of 850