



| | |
|---------------|-------------------|
| Product Code | 32124 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|---------|---------------|
| UK-REP | Not required |
| CH-REP | Not required |
| EAN/UPC | 0732094025163 |

Description

Welch Allyn bulbs:

- W.A. ref: 03100-U
- Compatible with following Welch Allyn product:
 - Otoscopes: 25020, 21700, 20200
 - Handle: 73500, holder: 28100
 - Illuminators: 27000, 27050, 26530, 41100, 43300
 - Old otoscopes: 25000, 25200, 20000