



Product Code	32152
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC 4045396219366

## Description

RCS-100 accessories: Endoscope adapter 4°

Riester code: 13277

## Technical Specifications

**Compatibility:** Compatible with flexible and rigid endoscopes with diameters from 2.7mm to 10mm (two versions planned, one optimized for optics between 2.7-4mm and one for 4-10mm) from a broad range of leading manufacturers based on C-mount standard (e.g. Wolf, Storz, Olympus, Scholly)

**Quick coupling:** Adapter can be attached to RCS-100 and a C-Mount endoscope quickly and it's ready to use