



Product Code	32153
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC	4045396213463
---------	---------------

Description

RCS-100 accessories: Endoscope adapter 9°

Riester code: 13278

Technical Specifications

Compatibility: Compatible with flexible and rigid endoscopes with diameters from 2.7mm to 10mm (two versions planned, one optimized for optics between 2.7-4mm and one for 4-10mm) from a broad range of leading manufacturers based on C-mount standard (e.g. Wolf, Storz, Olympus, Scholly)

Quick coupling: Adapter can be attached to RCS-100 and a C-Mount endoscope quickly and it's ready to use