



Product Code	32155
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC	4045396212053
---------	---------------

Description

RCS-100 accessories: General lens

Riester code: 13271

Technical Specifications

Focal length/width: 2.0
 Field angle: 78°
 Depth of field scope: 30mm ~ 4m
 Lighting source: natural light LED
 Magnification optical: 2x / digital: 4x
 LED colour temperance: 5500k
 Size - weight: 60.5x 60.5x19.0 mm - 66.4 g