



Product Code	32156
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2729185
CND	Z12040199
EAN/UPC	4045396212060

Description

RCS-100 accessories: Otoscope lens with pneumatic otoscopy port

Riester code: 13272

Technical Specifications

Focal length/width: 2.9

Field of view:

- 12mm Ø field of view at 25mm distance from viewing object (with 4mm specula size)
- 9mm Ø field of view at 15mm distance from viewing object (with 4mm specula size)
- 8mm Ø field of view at 10mm distance from viewing object (with 4mm specula size)

Max object distance: 15mm, at max object distance FOC Ø: 15mm

Depth of field scope: 10mm

Lighting source: natural light LED

Magnification: display: 6x

LED colour temperance: 4000 K

Size - weight: 73.5x40.0x40.0 mm - 96 g