



Product Code	32163
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2570722
CND	Z12021009
EAN/UPC	4719884250291
GMDN	12849

## Description

### ELITE PACK MS102

#### MD Scope with

- 1 camera probe 47 mm Ø 4.2 mm
- mini USB cable set
- 12 disposable specula
- internal Li-ion battery and charging cradle
- carrying case

#### Instant View

A clear and sharp image can be acquired and displayed instantly on a built-in screen. Just press a single button and wait for 3 seconds.

#### Real Time Interaction

Medical image can be captured and shared simultaneously to improve interaction between health professionals and patients.

#### Electronic medical record

A medical image is automatically generated to support most of the stand-alone health information system. Simply plug the cable into the video-out socket to retrieve the image.