

FLEXIBLE CAMERA PROBE 150 x diam 5.2 mm for MD Scope



Product Code	32168
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	Medical device
------	----------------

Class	I
-------	---

UK-REP	No
--------	----

CH-REP	No
--------	----

RDM (NSIS)	2804711
------------	---------

CND	Z12021009
-----	-----------

EAN/UPC	8023279321685
---------	---------------

GMDN	12849
------	-------

Description

MD SCOPE VIDEO OTOSCOPE ACCESSORIES AND SPARE PARTS:
Flexible camera probe 150 x Ø 5.2 mm