

LONG FLEXIBLE CAMERA PROBE 300 x diam 5.2 mm for MD Scope



Product Code	32169
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2804711
CND	Z12021009
EAN/UPC	8023279321692
GMDN	12849

Description

MD SCOPE VIDEO OTOSCOPE ACCESSORIES AND SPARE PARTS:
Long flexible camera probe 300 x Ø 5.2 mm