

MD SCOPE VIDEO OTOSCOPE MS102 - DELUXE ZIPPER PACK 2 PROBES



Product Code	32170
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	2570722
CND	Z12021009
EAN/UPC	4719884250383

Description

DELUXE PACK MS102 - 2 cameras
MD Scope supplied with 2 camera probes 47 (Ø 4.2 mm)/150 - cable set: RCA cable, charging cable - 12 disposable specula - internal Li-ion battery
Instant View: a clear and sharp image can be acquired and displayed instantly on a built-in screen. Just press a single button and wait for 3 seconds. Real Time Interaction: medical image can be captured and shared simultaneously to improve interaction between health professionals and patients. Electronic medical record: a medical image is automatically generated to support most of the stand-alone health information system. Simply plug the cable into the video-out socket to retrieve the image.

Technical Specifications

Image resolution: VGA (640 x 480 pixels)
LCD display: 2.8" panel
LEDs: 4 LEDs
Image capture: yes
Video recording: -
File storage: up to 1,000 photos
Video output: mini USB
File transfer/storage: mini USB
Time setting: -