

MD SCOPE VIDEO OTOSCOPE MS102 - DELUXE ZIPPER PACK 4 PROBES



Product Code	32171
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	2570722
CND	Z12021009
EAN/UPC	4719884250390

Description

DELUXE PACK MS102 - 4 cameras

MD Scope supplied with

- 4 camera probes 47 (Ø 4.2 mm)/75/150/300 mm
- cable set: RCA cable, charging cable
- 12 disposable specula
- internal Li-ion battery

Instant View: a clear and sharp image can be acquired and displayed instantly on a built-in screen. Just press a single button and wait for 3 seconds.

Real Time Interaction: medical image can be captured and shared simultaneously to improve interaction between health professionals and patients.

Electronic medical record: a medical image is automatically generated to support most of the stand-alone health information system. Simply plug the cable into the video-out socket to retrieve the image.

Technical Specifications

Image resolution: VGA (640 x 480 pixels)

LCD display: 2.8" panel

LEDs: 4 LEDs

Image capture: yes

Video recording: -

File storage: up to 1,000 photos

Video output: mini USB

File transfer/storage: mini USB

Time setting: -