

SILICONE SPECULA diam. 4.2 mm - reusable for 32166



Product Code	32173
Unit of sale	box of 24 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2804711
CND	Z12021009
EAN/UPC	8023279321739
GMDN	35348

Description

MD SCOPE VIDEO OTOSCOPE ACCESSORIES AND SPARE PARTS:
Reusable silicone specula Ø 4.2 mm for code 32166 - box of 24