

# MIC Wi-Fi & USB POLARIZED+UV+WHITE LED DERMATOSCOPE with software



Product Code	32177
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2184077
CND	Z12040118
EAN/UPC	8050713880132

## Description

WI-FI - USB DIGITAL POLARIZED+UV+WHITE LED DERMATOSCOPE with software

Dermatoscope with polarized filter is a very innovative handheld colour digital device with USB and WI-FI transmission. Suitable for skin and scalp. Small, easy to operate and portable.  
3 types of lighting: white, UV and white polarized LEDs.

Magnification:  
60x and 200x with 7 mm spacer ring or closed cap  
30x with transparent tube

MDR 2017/745/UE - Medical Device Class I  
Mic-Fi Derma software for Windows (GB, IT) and MAC (GB) to capture photos, compare images, history and reporting, data base for patients and selection of the point where analysis is carried out.

## Technical Specifications

Operating mode: continuous  
Frame rate: 15 frame/sec.  
Camera: Colour CMOS 1/4"  
Image format: bmp, jpg  
LEDs: 4 white, 4 UV and 4 polarized adjustable LEDs  
USB: USB 2.0, 2 mt  
Battery: Li-ion  
Working time: 2 hours  
Lens: high definition microscopy lens  
Power supply: DC 5.0 V/1A  
Consumption: Max 2.5 W  
Resolution: 1.3 Mpixel  
Size - weight: 36x142 mm - 88 g

## Standard accessories

USB cable  
Ruler  
1 Bag L + 1 bag S  
USB AC Adapter  
Contact closed cap  
Quick guide + manual (GB, FR, IT, ES, DE)  
USB foot pedal

