



Product Code	32179
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8050713880781

Description

DIGITAL WI-FI IRISCOPE with software and stand

Innovative handheld colour digital iriscopes with Wi-Fi wireless transmission. Small, easy to operate, light and portable.

2 adjustable white LEDs and 2 adjustable yellow LEDs.
It can be used with all the most popular smartphones, tablet, PC.

Iridology professional software (USB mode only) compatible with Windows 7, 8, 10, 11, available in GB, IT.

Technical Specifications

Magnification: 30x
 Camera: CMOS 1/2.7"
 Light source: 2 white leds and 2 yellow leds
 LED system: 3.2 V/20 ma
 Resolution: 5 Mpixel in USB mode
 Lens: high definition microscopy lens
 Frame rate: 15 frames/second (1280x1024), 30 frames/second (640x480) 320x240
 Power supply: Li-ion battery, DC 5.0V/1A
 Working time: continuous working, 2 hours fully charged
 Power consumption: 2.5 W
 Size - weight: Ø 36x142 mm - 88 g

Standard accessories

USB AC adapter
 USB cable
 Ring
 Chinrest
 Quick manual: GB, IT
 Software for Windows 7, 8, 10, 11 (GB, IT)