



Product Code	32180
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2183410
CND	Z12021009
EAN/UPC	8023279321807

## Description

Wi-Fi & USB VIDEOTOSCOPE with software

Innovative, handheld colour digital Otoscope with Wi-Fi wireless transmission. Small, easy to operate, light and portable. It can be used with all the most popular smartphones, tablets and even with PC.

## Technical Specifications

Operation system: Windows (7, 10, 11), Mac (OS X 10.7 and following), iOS, Android  
 Photo and video resolution: 1.3 MPixels - 1280x1024 (MJPG), 640x480 (MJPG), 320x240  
 Frame rate: 15 fps 1280x1024, 30 fps 640x480 320x240  
 Wi-Fi wireless transmission distance: up to 5 m  
 Power consumption: max 2.5 W  
 Li-ion battery: continuous working time - approx. 2 hours, full charging time - approx. 2 hours  
 Power source: DC 5.0 V / 1A  
 Size: Ø 36x142 mm  
 Weight: 88 g

Camera  
 Camera sensor: CMOS 1/4"  
 Optical features: Exposure: Auto, White balance: Auto, Gamma: Auto, Black reference: Auto

Otoscope  
 Magnification: 20 X  
 Focus range: 8-60 mm  
 Object distance: 1-160 mm  
 Light source: 4 adjustable white LED  
 Lens: high definition microscopy lens

Software interface: Mic-Fi Apps for Android and iOS system and Mic-Fi softwares for PC and MAC

## Standard accessories

9 disposable speculum (3x3.5, 3x4.3, 3x5.0 mm)  
 Plastic stand  
 USB AC adapter and cable  
 Manual (GB, FR, IT, ES, DE)  
 CD basic software Win (GB, FR, IT, ES, DE, UA, JP) Mac (GB)