



Product Code	32184
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279321845

Description

IRISCOPE Wi-Fi
Pour l'analyse de l'iris.

An innovative line of handheld colour digital microscopes: dermatoscope with polarized filter, capillaroscope and iridoscope with Wi-Fi wireless transmission. Small, easy to operate, light and portable. They can be used with all the most popular smartphones, tablet and PC. Basic software for acquisition of photos and videos. Google/Apple applications included.

Technical Specifications

IRISCOPE SPECIFICATIONS

Magnification: 30X

Light source: 2 adjustable white LED +2 adjustable yellow LED

Lens: High definition microscopy lens

Operation system Windows (7, 10, 11), OS X, iOS, Android (it can be used by Wi-Fi or USB)

Photo and video resolution 1.3 MPixels - 1280x1024 (MJPEG), 640x480 (MJPEG), 320x240

Frame rate 15 fps 1,280x1,024, 30 fps 640x480 320x240

Wi-Fi wireless transmission distance Up to 5 m

Power consumption Max 2.5 W

Li-ion battery Continuous working time - appx. 2 hours, full charging time - appx. 2 hours

Power source - Size - Weight DC 5.0V / 1A - Ø 36x142 mm - 88 g

CAMERA

Camera sensor CMOS 1/4"

Optical features Exposure: Auto, Gamma: Auto, White balance: Auto, Black reference: Auto

Standard accessories

USB AC adapter

USB cable

Quick manual

Manual (GB, FR, IT, ES, DE)

CD basic software Win (GB, FR, IT, ES, DE, UA, JP). Mac (GB)