



Product Code	32185
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	0806808145246

Description

MICFIEYE - Wi-Fi/USB CAMERA (C-mount)

MICFIEYE is a Wi-Fi camera with C-mount and allows to digitalize the image on tablets, smartphones and PC using dedicated Apps and Softwares.

MICFIEYE needs an adapter to be compatible with all optical devices such as microscopes, endoscopes, boroscopes, fiberscopes, C-mount lenses etc.

Supplied with usb cable (1.2 m), USB AC adapter, user manual (GB, IT, FR, ES, DE), software can be downloaded free of charge from web

Technical Specifications

Camera Type: Colour CMOS camera, 1.3 MP resolution

Sensor: ¼ Imax CMOS

Image pixel: 320x240, 640x480, 1280x1024

Lens standard: C-CS mount CCTV lens

Automatic optical features: Gamma correction, exposure, white balance, black reference, flicker avoidance, defect identification, colour saturation

Button snapshot: yes for Wi-Fi and USB mode

Dimension: Ø 36x75 mm

Wi-Fi transmission distance: =5 m

Video resolution: 1.3 MP, 320x240, 640x480, 1280x1024

Frequency range: 2,400-24,835 GHz

Li-Ion battery: Continuous working time - approx 2.5 hours

Power source: 5 V / 1A or USB

OS support: iOS / Android / Windows / Mac

Software features: Capture picture video recording, calibration and measurement.