

LITTMANN CLASSIC III - 5646 - copper - black



Product Code	32234
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

GTIN	00707387815233
------	----------------

Description

3M™ LITTMANN® CLASSIC III STETHOSCOPE SPECIAL EDITION

- 3M Code: 5646
- Colour: black - copper
- Lenght: 69 cm
- Head Ø: 43/33 mm

The 3M™ Littmann® Classic III™ Stethoscope offers high acoustic sensitivity for exceptional performance when doing general physical assessments. It features dual tunable diaphragms and updated design.

- versatile adult and pediatric two-sides chestpiece
- dual single-piece tunable diaphragms (pressure sensitive)
- small-side diaphragm prevents debris in open bell
- high acoustic sensitivity
- snap tight soft-sealing eartips
- next generation longer-life tubing
- small side converts to open bell